

# Референции

## Структура на OML\_021 съобщението

SEQ	OPT	RPT / #	GROUP	NAME
1	R	1		MSH Message Header
2	O	*		SFT Software Segment
3	O	*		NTE Notes and Comments
	O	1	Y	- PATIENT
4	R	1		<b>PID</b> Patient Identification
5	O	1		PD1 Patient Additional Demographic
6	O	*		NTE Notes and Comments
7	O	*		NK1 Next of Kin / Associated Parties
	O	1	Y	- PATIENT VISIT
8	R	1		<b>PV1</b> Patient Visit
9	O	1		PV2 Patient Visit - Additional Information
	O	*	Y	- INSURANCE
10	R	1		IN1 Insurance
11	O	1		IN2 Insurance Additional Information
12	O	1		IN3 Insurance Additional Information, Certification
13	O	1		GT1 Guarantor
14	O	*		AL1 Patient Allergy Information
	R	*	Y	- ORDER
15	R	1		ORC Common Order
	O	*	Y	- TIIMING
16	R	1		TQ1 Timing/Quantity
17	O	*		TQ2 Timing/Quantity Relationship
	O	1	Y	- OBSERVATION REQUEST
18	R	1		OBR Observation Request
19	O	1		TCD Test Code Detail
20	O	*		NTE Notes and Comments
21	O	1		CTD Contact Data
22	O	*		DG1 Diagnosis
	O	*	Y	- OBSERVATION
23	R	1		OBX Observation/Result
24	O	1		TCD Test Code Detail
25	O	*		Notes and Comments
	O	*	Y	- SPECIMEN
26	R	1		SPM Specimen
27	O	*		OBX Observation/Result
	O	*	Y	- CONTAINER
28	R	1		SAC Specimen Container detail
29	O	*		OBX Observation/Result
	O	*	Y	- PRIOR RESULT

SEQ	OPT	RPT / #	GROUP	NAME
	O	1	Y	- PATIENT PRIOR
30	R	1		PID Patient Identification
31	O	1		PD1 Patient Additional Demographic
	O	1	Y	- PATIENT VISIT PRIOR
32	R	1		PV1 Patient Visit
33	O	1		PV2 Patient Visit - Additional Information
34	O	*		AL1 Patient Allergy Information
	R	*	Y	- ORDER PRIOR
35	O	1		ORC Common Order
36	R	1		OBR Observation Request
37	O	*		NTE Notes and Comments
	O	*	Y	- TIMING PRIOR
38	R	1		TQ1 Timing/Quantity
39	O	*		TQ2 Timing/Quantity Relationship
	R	*	Y	- OBSERVATION PRIOR
40	R	1		OBX Observation/Result
41	O	*		NTE Notes and Comments
42	O	*		FT1 Financial Transaction
43	O	*		CTI Clinical Trial Identification
44	O	1		BLG Billing

From:

<https://hl7.skyware-group.com/> - **iLab HL7 Wiki**

Permanent link:

<https://hl7.skyware-group.com/doku.php?id=full-ref&rev=1506087570>

Last update: **2017/09/22 16:39**

