

Резултат (ORU_R01)

Съобщението ORU_R01 (Unsolicited transmission of an observation message) се изпраща от ЛИС към БИС. Изпраща се след валидация на резултатите във визита в ЛИС.

Структура на съобщението

По долу е дадена **редуцирана** структура на ORU_R01 съобщението, която съдържа само релевантните към настоящата интеграция групи и сегменти. Наличие на други валидни групи и сегменти няма да попречи на обработката, те ще бъдат игнорирани.

Сегмент / Група	Кол. по стандарт	Кол. за iLab
MSH Message Header	1	1
SFT Software Segment	0..*	0..*
PATIENT RESULT начало	1..*	1
PATIENT начало	0..1	1
PID Patient Identification	1	1
NTE - Notes and Comments	0..*	0..1
PATIENT_VISIT начало	0..1	1
PV1 Patient Visit	1	1
PATIENT_VISIT край		
PATIENT край		
ORDER OBSERVATION начало	1..*	1..*
ORC - Common Order	0..1	1
OBR - Observation Request	1	1
NTE - Notes and Comments	0..*	0..1
OBSERVATION начало	0..*	0..*
OBX - Observation/Result	1	1
NTE - Notes and Comments	0..*	0..1
OBSERVATION край		
SPECIMEN начало	0..*	0..*
SPM - Specimen	1	1
SPECIMEN край		
ORDER OBSERVATION край		
PATIENT RESULT край		

Пример

```
MSH|^~\&|iLab|Лаборатория|iHealth|МБАЛ АД|20171104230436||ORU^R01^ORU_R01|2a3fdb1de|P|2.5.1|
SFT|SKYWARE Group|5.0.0.41525|iLabCourier|5.0.0.41525|
PID|1||9706166900^^^GRAO^NI~16403^^^Laboratory^MR~9598/2017^^^Hospital^MR~5221^^^LIS^X
X||Бенков^Димитър^Иванов||19970616|M|
PV1|1|N|Урологични легла^^^01.03|||||||16403^^^Laboratory^MR|
```

ORC|RE|IH20171028103353461|16403^iLab|||^201710281041|||1999900024^Русинов^Атанас^Марин
ов^^Д-р^^BLS^^^DN|

OBX|1|IH20171028103353461|16403^iLab|57021-8^Пълна Кръвна
Картина^LN^1-2^HCPT||201802121826

OBX|1|NM|6690-2^Левкоцити^LN^0-3^HCPT||5.17|10\9/L|3.50 -
10.50|||C||201710281036||2012312345^Димитрова^Анета^^^Д-р^^BLS^^^DN|

OBX|12|NM|789-8^Еритроцити^LN^0-4^HCPT||3.00|10\12/L|4.20 -
5.40|L||C||201710281036||2012312345^Димитрова^Анета^^^Д-р^^BLS^^^DN|

OBX|13|NM|718-7^Хемоглобин^LN^0-1^HCPT||91|g/L|125 -
180|LL||C||201710281036||2012312345^Димитрова^Анета^^^Д-р^^BLS^^^DN|

OBX|14|NM|4544-3^Хематокрит^LN^0-2^HCPT||0.273|L/L|0.380 -
0.540|L||C||201710281036||2012312345^Димитрова^Анета^^^Д-р^^BLS^^^DN|

OBX|15|NM|787-2^MCV^LN^0-5^HCPT||91.0|fL|82.0 -
92.0|||C||201710281036||2012312345^Димитрова^Анета^^^Д-р^^BLS^^^DN|

OBX|16|NM|785-6^MCH^LN^0-6^HCPT||30.5|pg|28.0 -
33.0|||C||201710281036||2012312345^Димитрова^Анета^^^Д-р^^BLS^^^DN|

OBX|17|NM|786-4^MCHC^LN^0-7^HCPT||335|g/L|320 -
360|||C||201710281036||2012312345^Димитрова^Анета^^^Д-р^^BLS^^^DN|

OBX|18|NM|777-3^Тромбоцити^LN^0-10^HCPT||240|10\9 /L|140 -
440|||C||201710281036||2012312345^Димитрова^Анета^^^Д-р^^BLS^^^DN|

SPM|1||BLD^кръв^HL70487|||201710281036|

ORC|RE|IH20171028103353461|16403^iLab|||^201710281041|||1999900024^Русинов^Атанас^Марин
ов^^Д-р^^BLS^^^DN|

OBX|2|IH20171028103353461|16403^iLab|4537-7^СУЕ^LN^0-20^HCPT||201802121826

OBX|1|NM|4537-7^СУЕ^LN^0-20^HCPT||75|mm/h|3 -
46|H||C||201710281036||2012312345^Димитрова^Анета^^^Д-р^^BLS^^^DN|

SPM|1||BLD^кръв^HL70487|||201710281036|

ORC|RE|IH20171028103353461|16403^iLab|||^201710281041|||1999900024^Русинов^Атанас^Марин
ов^^Д-р^^BLS^^^DN|

OBX|3|IH20171028103353461|16403^iLab|14749-6^Глюкоза^LN^0-9^HCPT||201802121826

OBX|1|NM|14749-6^Глюкоза^LN^0-9^HCPT||9.89|mmol/l|3.30 -
6.40|H||C||201710281036||2012312345^Димитрова^Анета^^^Д-р^^BLS^^^DN|

SPM|1|^AITG||SER^серум^HL70487|||201710281036|

From:
<https://hl7.skyware-group.com/> - **iLab HL7 Wiki**

Permanent link:
<https://hl7.skyware-group.com/doku.php?id=oru-r01>

Last update: **2018/09/26 19:06**

